




POST AND STAND
NOT SHOW IN THIS VIEW



POST AND STAND
BASE DETAIL

NOTE:

$\triangle 1$ DETECTOR SURFACE IS 5.5mm (0.22in) BELOW FRONT FACE OF APERTURE PLATE.

		5100 Patrick Henry Drive Santa Clara, CA 95054 www.coherent.com	
Title: TOP ASSEMBLY, PowerMax USB-PM2X		PART NUMBER: 1257617	
 THIRD ANGLE PROJECTION		PROPRIETARY NOTICE: THIS DOCUMENT IS THE PROPERTY OF COHERENT, INC. AND CONTAINS PROPRIETARY INFORMATION AND IS DELIVERED ON THE EXPRESS CONDITION THAT IT AND THE INFORMATION CONTAINED THEREIN ARE NOT TO BE USED, DISCLOSED OR REPRODUCED IN WHOLE OR IN PART FOR ANY PURPOSE (INCLUDING MANUFACTURE FOR ANYONE OTHER THAN COHERENT) WITHOUT THE EXPRESS WRITTEN CONSENT OF COHERENT.	
DRAWING NOT TO SCALE		Date: 2-15-19	Rev 01