





COOLANT FITTINGS INCLUDED
1/8" NPT PORTS IF FITTINGS REMOVED.




2X M4X0.7 - 6H ∇ .32in [8.1mm]
MOUNTING THREADS FOR
INCLUDED MOUNTING PLATE
SEE SHEET 3 OF 3



\square ϕ .17in [4.3mm] ∇ .13in [3mm]
2X M4X0.7 - 6H ∇ .44in [11mm]
MOUNTING THREADS FOR
INCLUDED MOUNTING PLATE
SEE SHEET 3 OF 3




PowerMax-Pro 1 kW Average Models				
PART NUMBER	DESCRIPTION	MAXIMUM AVERAGE POWER	MAXIMUM PEAK POWER	WINDOW INCLUDED
1324794	PMP kW - 1kW 810NM	1kW	1kW	NO
1324795	PMP kW - 1kW 10.6 μ m	1kW	1kW	NO
1324796	PMP kW - 1kW 1.1 μ m	1kW	1kW	NO
1324797	PMP kW - 1kW 810NM, 1.1 μ m 10.6 μ m	1kW	1kW	NO
1325221	PMP kW - 1kW 1.1 μ m	1kW	1kW	YES



5100 Patrick Henry Drive
Santa Clara, CA 95054
www.coherent.com

Title:
**POWERMAX-PRO 1 KILOWATT HD
SENSOR**

THIRD ANGLE
PROJECTION



PART NUMBER:

PROPRIETARY NOTICE:
THIS DOCUMENT IS THE PROPERTY OF COHERENT, INC. AND CONTAINS PROPRIETARY INFORMATION AND IS
DELIVERED ON THE EXPRESS CONDITION THAT IT AND THE INFORMATION CONTAINED THEREIN ARE NOT TO BE
USED, DISCLOSED OR REPRODUCED IN WHOLE OR IN PART FOR ANY PURPOSE (INCLUDING MANUFACTURE FOR
ANYONE OTHER THAN COHERENT) WITHOUT THE EXPRESS WRITTEN CONSENT OF COHERENT.




DRAWING NOT TO SCALE

Date:

Rev 01




2X M4X0.7 - 6H ∇ .32in [8.1mm]
MOUNTING THREADS FOR
INCLUDED MOUNTING PLATE
SEE SHEET 3 OF 3



5100 Patrick Henry Drive
Santa Clara, CA 95054
www.coherent.com

Title:
POWERMAX-PRO 1 KILOWATT HD
SENSOR

PART NUMBER:








PROPRIETARY NOTICE:
THIS DOCUMENT IS THE PROPERTY OF COHERENT, INC. AND CONTAINS PROPRIETARY INFORMATION AND IS
DELIVERED ON THE EXPRESS CONDITION THAT IT AND THE INFORMATION CONTAINED THEREIN ARE NOT TO BE
USED, DISCLOSED OR REPRODUCED IN WHOLE OR IN PART FOR ANY PURPOSE (INCLUDING MANUFACTURE FOR
ANYONE OTHER THAN COHERENT) WITHOUT THE EXPRESS WRITTEN CONSENT OF COHERENT.

DRAWING NOT TO SCALE

Date:

Rev 01



 COHERENT 5100 Patrick Henry Drive Santa Clara, CA 95054 www.coherent.com		PART NUMBER:	
Title: POWERMAX-PRO 1 KILOWATT HD SENSOR		Rev 01	
 THIRD ANGLE PROJECTION		DRAWING NOT TO SCALE	
Date:		Rev 01	