

NOTE:

1 DETECTOR SURFACE IS 18.4mm (0.73in) BELOW FRONT FACE OF APERTURE PLATE

Part Number	Description
1232163	PowerMax USB - PM1K Power Sensor
1276824	PowerMax USB - PM3K Power Sensor
1215999	PowerMax USB - PM5K Power Sensor




**COHERENT**

5100 Patrick Henry Drive  
Santa Clara, CA 95054  
www.coherent.com

Title:  
PowerMax USB - PM1K/PM3K/PM5K Power  
Sensor

PART NUMBER:  
-



THIRD ANGLE  
PROJECTION

PROPRIETARY NOTICE:  
THIS DOCUMENT IS THE PROPERTY OF COHERENT, INC. AND CONTAINS PROPRIETARY INFORMATION AND IS  
DELIVERED ON THE EXPRESS CONDITION THAT IT AND THE INFORMATION CONTAINED THEREIN ARE NOT TO BE  
USED, DISCLOSED OR REPRODUCED IN WHOLE OR IN PART FOR ANY PURPOSE (INCLUDING MANUFACTURE FOR  
ANYONE OTHER THAN COHERENT) WITHOUT THE EXPRESS WRITTEN CONSENT OF COHERENT.

DRAWING NOT TO SCALE Date: 6-23-15 Rev 01